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Thank you for purchasing a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long-term and flawless operation.

SUMMARY

SUMMARY	42
STANDARD ITEMS AND OPTIONS	44
GENERAL	45
SAFETY INFORMATION	45
A. SET-UP	46
1. GENERAL	46
2. INSTALL PERISTALTIC PUMP HEADS ON PUMP DRIVES:	46
2.1. SP QUICK.....	46
2.2. SP QUICK D (2 EA. SP QUICK PUMP HEADS ARRANGED IN SERIES)	47
2.3. SP STANDARD, SP VARIO	47
3. INSTALL HOSE - PERISTALTIC PUMPS.....	48
3.1. GENERAL	48
3.2. HOW TO CHOOSE THE HOSE MATERIAL.....	48
3.3. SP QUICK.....	48
3.4. SP QUICK D.....	49
3.5. SP STANDARD	49
3.6. SP VARIO	49
3.7. USE OF SMALL ID HOSES (SIZE 0.8 AND 1.7 MM).....	50
4. INSTALL MULTI-CHANNEL PUMP HEAD ON PUMP DRIVE	50
4.1. INSTALL MULTI-CHANNEL ADAPTER ON PUMP DRIVE	50
4.2. INSTALL MULTI-CHANNEL PUMP HEAD ON PUMP DRIVE	51
5. INSTALL HOSE AND CASSETTE	52
5.1. CASSETTE "SMALL"	52
5.1.1. HOW TO OPTIMIZE CONTACT PRESSURE SETTING.....	53
5.1.1.1. For maximum of hose life	53
5.1.2. PUMP OPERATION WITH PARTIAL USE OF A MULTI-CHANNEL PUMP HEAD	53
5.1.3. UNLOAD HOSE CASSETTE:	54
5.1.4. STOPPER HOSE EXTENSION:	54
5.1.4.1. Size 0.2; 0.5; 0.9 and 1.4 hoses	54
5.1.4.2. Size 2.8 hose	55
5.2. CASSETTE MEDIUM;	55
5.3. CASSETTE LARGE.....	55
B. OPERATION	57
6. GENERAL	57
7. PD 5001, PD 5006	57
8. PD 5101; PD 5106	57
8.1. INTERFACE	59
9. PD 5201, PD 5206	60
9.1. STANDARD FUNCTIONS	60
9.2. CHARGING / DISCHARGING OF HOSE	61
9.3. EXTRA FUNCTIONS	61
9.3.1. HOW TO PROGRAM PARTICULAR PUMP HEADS.....	62
9.3.2. SELECT DISPLAY MODE	62
9.3.3. DISPLAY FEED RATE IN ML/MIN.....	62
9.3.4. CALIBRATE FEED RATE.....	63
9.3.5. METERING BY VOLUME	63
9.3.6. CALIBRATE METERED VOLUME.....	64

9.3.7.	INTERVAL METERING.....	65
9.3.8.	ADJUST BACKSTROKE FOR METERING BY VOLUME MODE.....	65
9.4.	INTERFACES	65
9.4.1.	ANALOG INTERFACE	65
9.4.2.	RS 232 INTERFACE	66
9.5.	RECOMMENDED SIZES OF HOSE	67
C. RECOMMENDATIONS ABOUT HOSES		68
10.	HOSE PROPERTIES	68
11.	CHEMICAL STABILITY OF HOSES.....	70
12.	SERVICE LIFE OF HOSES	72
13.	FEED RATES - LIST OF HOSES FOR PERISTALTIC PUMPS	73
14.	FEED RATES - LIST OF HOSES FOR MULTI-CHANNEL PUMPS	74
CLEANING AND SERVICING.....		75
UNINSTALL, FORWARD & STORE.....		75
DISPOSAL.....		75
TROUBLESHOOTING		76
SPECIFICATIONS.....		77
15.	PERISTALTIC PUMP DRIVES	77
16.	FEED RATE GRAPHS FOR PERISTALTIC PUMPS	77
17.	FEED RATE GRAPHS FOR MULTI-CHANNEL PUMPS	78
WARRANTY, LIABILITY & COPYRIGHT		78
FAQ / REPAIR WORK		78
CE-DECLARATION OF CONFORMITY		79

Important information

Advice about power cord / mains supply

Caution: mandatory action

Caution: fire- and explosion hazard

Advice about maintenance / repair

STANDARD ITEMS AND OPTIONS

	item	quantity	P/N 230V / 50/60Hz	P/N 115V / 50/60Hz
	PD 5001 pump drive	1	523-50010-00	523-50010-01
or	PD 5001 pump drive includes multi-channel adapter	1	523-50013-00	523-50013-01
or	PD 5006 pump drive	1	523-50060-00	523-50060-01
or	PD 5101 pump drive	1	523-51010-00	523-51010-01
or	PD 5101 pump drive includes multi-channel adapter	1	523-51013-00	523-51013-01
or	PD 5106 pump drive	1	523-51060-00	523-51060-01
or	PD 5201 pump drive	1	523-52010-00	523-52010-01
or	PD 5201 pump drive includes multi-channel adapter	1	523-52013-00	523-52013-01
or	PD 5206 pump drive	1	523-52060-00	523-52060-01
or	PD 5001 with SP quick 1,6	1	523-50019-00	523-50019-01
or	PD 5006 with SP quick 1,6	1	523-50069-00	523-50069-01
or	PD 5006 with SP standard 2,5	1	523-50068-00	523-50068-01
or	PD 5101 with SP quick 1,6	1	523-51019-00	523-51019-01
or	PD 5106 with SP quick 1,6	1	523-51069-00	523-51069-01
or	PD 5106 with SP standard 2,5	1	523-51068-00	523-51068-01
or	PD 5201 with SP quick 1,6	1	523-52019-00	523-52019-01
or	PD 5206 with SP quick 1,6	1	523-52069-00	523-52069-01
or	PD 5206 with SP standard 2,5	1	523-52068-00	523-52068-01
	PD 5201 with C4 und 4 x Cassette small	1	523-52017-00	523-52017-01
	Instruction manual	1	01-005-002-75	01-005-002-75
	power cord (for PD 5001, PD 5006 only)	1	14-007-003-81	14-007-003-89

Accessories

	item	P/N
	peristaltic pump head SP quick (SWS 1.6 mm)	527-11100-00
	peristaltic pump head SP quick (SWS 2.5 mm)	527-11300-00
	peristaltic pump head SP quick d (SWS 1.6 mm)	527-11120-00
	peristaltic pump head SP quick d (SWS 2.5 mm)	527-11320-00
	peristaltic pump head SP standard (SWS 1.6 mm)	523-43010-00
	peristaltic pump head SP standard (SWS 2.5 mm)	523-43030-00
	peristaltic pump head SP vario	523-45110-00
	C 4 multi-channel pump head	524-80420-00
	C 8 multi-channel pump head	524-40810-00
	C 12 multi-channel pump head	524-81220-00
	cassette, small (SWS 0.9 mm)	524-90022-00
	cassette, medium (SWS 1.6 mm)	524-90021-00
	cassette, large (SWS 1.6 mm)	524-90010-00
	pedal switch	526-14100-00
	adapter for multi-channel pump heads	526-16000-00
	hose coupling 0.2-2,8 mm inner diameter hose	526-22000-00
	hoses as required, refer to Chapters 13 and 14	

GENERAL

-  Unpack your item carefully.
Inspect for damage and report such damage or missing parts to your supplier right away.
-  Read your Instruction Manual carefully. Take time to save time while working with your product. Make sure that every user has read and understood the Instruction Manual.
-  Please store the Instruction Manual in a place easily accessible to every user.

IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !

-  A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the products.
For the Continental US they feature a US-standard plug (NEMA Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).
-  For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one mating your local system.
-  As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !

SAFETY INFORMATION

Please comply with all safety and accident prevention regulations as in force for laboratory work !

Comply with all safety instructions in force for pumps !

Use extra care while working with flammable substances; refer to safety data sheets. Store safety data sheets in a place easily accessible to all operators.

Use eye protection and adequate clothing.

Use extra care while working with flammable substances. Refer to safety data sheets.

Pump drives do not feature explosion protection!

Caution ! Hoses are wear parts. Inspect hoses in regular intervals and replace on condition.

Before connecting your item with the local power supply, make sure the item has been designed for your local voltage; refer to data plate on item.

E

Turn power switch to OFF when item is out of use, or before disconnecting from mains.

Repair work is limited to skilled personnel authorized by Heidolph.

Your item needs a solid stand.

Before turning pump drives ON, make sure your hands are clear of the squeeze rollers; contusion hazard !

A. SET-UP

1. General (fig. 21)

- The pump drive is located on a smooth horizontal plane.
- When using more pump drives at the same time (2 ea. max.), they may be stacked to save space.

fig. 21

2. Install peristaltic pump heads on pump drives:

2.1. SP quick (figs. 22,23)

- Press SP quick pump head (3) against PD.... pump drive front panel so that thread (4) and coupling (5) are in line with threaded bolt and coupling on pump head.
- Protruding parts of the coupling crown need to slip into cut-outs provided in pump drive coupling stud. To this end tilt pump head around rotor center or move rotor with forefinger, hose bed being open (refer to Chapter 3.3).


fig. 22

- When coupling is seated properly, secure pump head (3) with 2 ea. knurled-head bolts (6) by hand; if required, tighten with screwdriver or coin.


fig. 23

2.2. SP quick d (2 ea. SP quick pump heads arranged in series) (figs. 24,23,22)

Limited to PD 5001, PD 5101 and PD 5201 pump drives.

- SP quick pump head with one more SP quick pump head allows to run 2 ea. feeder channels the same time.
- Unscrew 2 ea. knurled head bolts (6 fig. 23) from SP quick pump head (3) with about 10 revolutions CCW. Replace with long knurled head bolts (7).
- Beware short knurled head bolts (6 fig. 23) for future use.
- Fit SP quick d pump head (8), with spacer (42) over threaded end of long knurled head bolts (7), seat long knurled head bolts and secure by turning clockwise; press two ea. pump heads (3 & 8) together. SP quick (3) properly needs to be caught by the SP quick d (8) counterpiece. As required, rotate rotor by hand (to this end open hose bed, refer to Chapter 3.3) Secure pump pack on pump drive as described in Chapter 2.1.


fig. 24


fig. 25

2.3. SP Standard, SP vario, (fig 25)

- Installation of SP Standard and SP vario pump heads needs flange (9) from your hardware kit.
- Uninstall flange (9) from pump head (3) by loosening setscrew (10).
- Secure flange on front panel with countersunk screws (11). Use Phillips screwdriver.
- Introduce SP Standard or SP vario pump head into flange (coupling pointing forward).
- Make sure that coupling fits counterpiece (as required, rotate rotor a little bit; remove cover to get access (also refer to Chapter 3.5)). Check horizontal location of pump head (3); hose outlet pointing to the left; secure in place with setscrew (10).

3. Install hose - peristaltic pumps

3.1. General

To avoid damage to the hose, make sure you are going to install the right wall thickness. If not, pump may fail to feed. To determine the right wall thickness, make reference to table below:

wall thickness (mm)	1.6	2.5
SP quick	527-11100-00	527-11300-00
SP quick d	527-11120-00	527-11320-00
SP Standard	523-43010-00	523-43030-00
SP vario	523-45110-00	523-45110-00

3.2. How to choose the hose material

For more information about hoses, their properties and chemical stability make reference to Chapter C.

3.3. SP quick (fig. 26)

- The SP quick features snap-on mechanism to speed up installation of hoses.
- Fold lever (18) to the left to open hose bed.
- Rotate two clampscrews (19) CCW till stop.
- Seat hose (17) in sickle-shaped cut-out in pump head.
- Fold lever (18) to the right to close hose bed.
- Do a short trial run to check proper clamping of hose.
- Make shure that hose does not migrate during rotation.
- If necessary clamping can be increased by turning clampscrews (19) clockwise when hose bed is open. Close hose bed and check adjustment with trial run.
- Repeat above setting and testing work as required.


fig. 26

Carefully loosen item; remember that clamping power has not been adjusted before and excessive clamping power may squeeze-off the hose.

Check for hose clamps properly seated. Hose (17) should not creep while pumping; check hose clamping in a brief test run.

3.4. SP quick d (fig. 27)

Use is limited to PD 5001, PD 5101 and PD 5201 pump drives.

- If two pump heads are set-up in series, the hose will be installed in the forward pump head (3), as described in Chapter 3.3 SP quick. The hose itself will be seated in the rear pump head's hose bed (8) through a wide gap between both pump heads. Continue as described for SP quick in Chapter 3.3.


fig. 27

3.5. SP Standard (fig. 28)

- Loosen knurled nuts (22 & 23) to remove clamping lever (20) and cover (21). Seat hose (17) in hose bed (180 degrees loop), rotate rotor (24) by hand until hose (17) is pulled into the gap between the rollers (25) and hose bed (26).
- Check hose being squeezed properly by rotating the rotor (24) through several revolutions.
- Install clamp lever (20) and cover (21); hold both hose ends tight while tightening knurled nut (23).
- Holding hose ends tight over the rollers (25) will make sure both rollers (25) are centered over the hose (17) which is an essential requirement for proper feeding action and standstill sealing.
- Do a short test run (function test).


fig. 28

3.6. SP vario (fig. 29)

- The SP vario pump head accepts hoses of different wall thicknesses. Adjustment uses knurled knob (27).
- Set clockwise rotation with selector (31), refer to figs. 47, 49 and 51).
- Uninstall clamp lever (20), cover (21) by loosening knurled nuts (22 & 23).
- Set clearance to max. hose thickness with knurled knob (27); rotate in "+" direction.
- Introduce a 180 degrees hose (17) loop in hose bed, rotate rotor (24) by hand until hose (17) is pulled into the gap between rollers (25) and hose bed (26).
- Run rotor (24) through several revolutions until hose is squeezed over its entire width.
- Energize pump on condition.

Select "clockwise" rotation before energizing pump !

CAUTION ! Contusion hazard

- Hose thickness is adjusted with knurled knob (27). Turn to "+" increase clearance (thick wall hose) or "-" to decrease clearance (thin wall hose).
- Best setting is reached, when hose is squeezed over about 80 % of its wall thickness. At this point, the fluid filled hose is squeezed by the rollers as to allow feeding; the same time no fluid may drain from the hose at standstill.
- Install cover (21) and clamp lever (20); hold both hose ends tight while tightening knurled nut (23).
- For easy adjustment of different sizes of hose wall thicknesses use enclosed auxillary tool for radial hole on knurled screw (27). After adjustment tighten screw (44) to secure knurled screw.
- Pulling hose ends over the rollers (25) will make sure both rollers (25) are centered over the hose (17) which is an essential requirement for proper feeding action and standstill sealing.
- Do a short test run to verify setting. Repeat setting and test run on condition.
- Before energizing the pump, close cover (21); Contusion HAZARD !
- Depress lever (43) to uninstall rotor by pulling; refer to arrow mark for this action.


fig. 29


fig. 30

3.7. Use of small inner diameter hose (size 0.8 and 1.7 mm) (fig. 30).

- Short lengths of 0.8 and 1.7 hoses are supplied with the SP vario and SP Standard pump heads; they are used to improve clamping action.
- To this end slip these short hoses (28) over your pump hose (17) and install hose (17) with short hoses (28) in pump head so that these short hoses are located under the clamp. To go ahead from this point, refer to "install hose" Chapter.


fig. 31

Normally, the multi-channel adapter comes installed on the pump drive. If a multi-channel

adapter had been procured as an option for an existing pump drive, install item as stipulated below:

- Plug screws (50) from hardware kit in holes provided.
- Align multi-channel adapter coupling with counterpiece (51) on pump drive.
- Fit multi-channel adapter in coupling, align screws (50) with threaded holes (52) and pull both screws (50) tight with a screwdriver.

4.2. Install multi-channel pump head on pump drive (fig. 32)

Tooling required: long-reach screwdriver (from hardware kit).

- To install pump head on pump drive, align follower pin (54) with slot in multi-channel adapter (53, fig. 31).
- Secure pump head on pump drive with 4 ea. screws (55, from hardware kit); use special screwdriver from hardware kit.
- We recommend to support the pump head for this action; use adjustable foot (56) for leveling. Pull lock-nut (57) tight against the pump head's underside.
- CAUTION: the multi-channel pump head's drive C 8 follower pin features a plastic sleeve (58) over the follower pin (54). Check for proper seat of parts. The plastic sleeve (58) is used for cushion between both parts of the coupling, limiting noise emission to a minimum. This sleeve (58) needs to be inspected after every 500 hours of operation and replaced on condition. Sleeve P/N: 23-09-01-01-88-0


fig. 32


fig. 33

CAUTION: absence of SLEEVE (58) will harm the pump head.

5. Install hose and cassette

CAUTION: use only those combinations of multi-channel pump heads, cassettes, and hoses specified in reference table (Chapter 15). Deviation will increase wear or cause malfunction.

5.1. Cassette "small"

The cassette models "small" have been designed exclusively for C4 and C12 multi-channel pump heads.

The cassette model "small" accepts 3 only stopper hoses ! Use of hose in absence of a stopper will cause malfunction and voids your warranty !

Disconnect pump drive before installing or uninstalling cassettes, to avoid contusion hazard.

- The hose cassette system installs and uninstalls at ease (snap-on); no extra tooling required for this work. Since hose cassettes may work in both directions, no particular orientation of the cassette is required. Despite this fact, we recommend to arrange cassettes alternatingly, for ease of access to the setting knobs of every single cassette in a stack.
- Turn contact pressure setting knob (59) (fig. 34) to max. clearance (counter-clockwise); You do not need to do this work if you haven't touched the factory setting before, or contact pressure had been set for actual pumping work before.
- Load approved 3-stopper-hose in cassette (refer to fig. 35). Slip stopper on both ends of the cassette to the hose holder (refer to fig. 36).
- Hook hose cassette over one of the carrier rods, while holding hose between cassette feet and pressing its against the hose holder (refer to figs. 36 & 37).


fig. 34


fig. 35




fig. 36

- Fold hose over rotor, depress latch (refer to pig. 37) until pawl snaps in (refer to fig. 38).

Note: as required, depress hose cassette at place marked in fig. 38, to ease snap action.

- Adjust hose contact pressure. Turn knob (59) until inside edge of wedge is centered between #1 and #2 marks on the data plate graduation (refer to fig. 34).

5.1.1. How to optimize contact pressure setting

For certain applications, contact pressure needs re-adjustment . For example, varying feed rate of a particular hose, minimizing feed rate fluctuations caused by changes in system pressure, or to simply extension.

5.1.1.1. For maximum of hose life

With the pump running, decrease hose contact pressure adjusted before. To this end turn regulating knob (59, fig. 34) counterclockwise to displace wedges to a higher scale reading (towards 5). Go ahead turning knob (59, fig. 34) until feed rate or pressure are too low. From this point go ahead turning knob clockwise until feed rate and pressure are ok.

5.1.2. Pump operation with partial use of a multi-channel pump head (fig. 39)

Pump may be run with all or part of the multi-channel pump head's ports connected. Both metal rods of the pump feature an elastic packing (60) on both ends that secures cassettes on the pump. Secure multi-channel pump heads with part of the ports loaded with cassettes by slipping packings (60) over the rods towards the front end of the outboard cassette (refer to e.g. position A). On


fig. 39

E

pumps completely loaded with cassettes, make sure to slip the packings (60) against the front end of the pump as far as possible (refer to position B).

5.1.3. Unload hose cassette:

- Turn drive to OFF.
- Depress cassette's locking pawl (refer to fig. 38) and unhook from rod.
- Tilt cassette up and unhook from rod. Unload cassette by lifting.

5.1.4. Stopper hose extension:

Optional PTFE hose couplings are available for extension. Just plug together with hoses from the reel.


fig. 40

5.1.4.1. Size 0.2; 0.5; 0.9 and 1.4 hoses

Use packings from your hardware kit. Plug hose through coupling nut (refer to fig. 41); install packings and pull coupling nut hand-tight.


fig.41

5.1.4.2. Size 2.8 hose


fig.42

Use cone from your hardware kit. Refer to fig. 42 and plug hose through coupling nut, after having trimmed (pointed) hose end. After installation cut pointed end apart. Install cone as shown, pull coupling nut hand-tight.

5.2. Cassette medium;

Medium-model cassettes exclusively have been designed for C8 multi-channel pump heads.

Installation procedure is the same as described for the small-model cassette (refer to Chapter 5.1 ff); these cassettes however do not feature stoppers, they use hoses from the reel only (refer to Chapter 15, reference table).

When cassette is interlocked with multi-channel pump head exert some pull force on both ends of the hose to prevent hose from looping.

To adjust roller contact pressure, turn setting knob until inboard wedge reads between 3 and 4 on the scale; refer to Chapter 5.1.1 to optimize roller contact pressure.

5.3. Cassette large

Large-model cassettes exclusively have been designed for C8 multi-channel pump heads.

Installation procedure is the same as described for the small-model cassette (refer to Chapter 5.1 ff); these cassettes however do not feature stoppers, they use hoses from the reel only (refer to Chapter 15, reference table).

The large-model cassette features adjustable hose retainers on in- and outlet ends. These retainers adjust hose clamping action for different sizes of hoses. To adjust setting knob (61) make reference to fig. 43 and the reference table below:

size of hose	set to
1,7	14
3,1	16
4,8	25
6,3	17


fig. 43

E

When cassette is interlocked with multi-channel pump head exert some pull force on both ends of the hose to prevent hose from looping.

To adjust roller contact pressure, turn setting knob until inboard wedge reads between 3 and 4 on the scale; refer to Chapter 5.1.1 to optimize roller contact pressure.

B. OPERATION

6. General

Before running the item, read all of the safety and general information !

Before connecting your item with your local power supply, make sure item has been designed and wired for your voltage and frequency. Use your national standard plugs. For more information refer to Chapter "General Information".

- Connect with mains.

7. PD 5001, PD 5006 (fig. 44)

- Power switch (29) on control panel starts pump rotation.
- Select speed with selector knob (30).
- Feed rate is proportional with pump speed. Go by rpm scale on the item. For more information about exact speed refer to reference tables in Chapters 17 and 18.
- Select sense of rotation with switch (31); (go by arrow marks).

fig. 44

8. PD 5101; PD 5106 (figs 45, 46)

- PD 5101 and PD 5106 pump drives are splash proof (IP 55).

Service work is limited to expert personnel authorized to do so by Heidolph.

Keep interface connector (37) closed all the time by either an interface plug or a blind plug, not to void splash protection.

E

- Item's master switch (29) is located on the rear panel. LED (40) for sense of rotation in ON.
- Set speed with speed selector (30) on the control panel.
- Feed rate is proportional with pump speed. Go by rpm scale on the item. For more information about exact speed refer to "Specifications" section.

fig. 45


Bild 46

- 30: Rotation speed
- 31: Direction of rotation
- 32: Start/Stop
- 33: Maximal conveyance

- Select sense of rotation with switch (31); (go by arrow marks). Sense of rotation may be changed on standstill only.
- ⚠ Caution!** Hit Stop key (32) before changing sense of rotation, select new sense of rotation and restart.
- Sense of rotation selected is confirmed by LEDs (40).
 - Start pump with „Start Stop“ key (32). Pump status is confirmed by an LED (44).
 - To stop pump rotation, hit „Start Stop“ key (32) once more.
 - When hitting „max.“ key (33), pump will run at max. speed as long as this key is held depressed; this feature eases filling or draining hoses without any change to basic settings.

8.1. Interface (fig. 47)

- This built-in analog interface allows external pump control.

Control voltage 0 - 10 volts equals a rotational speed of 0 to max.

control voltage	PD 5101	PD 5106
0 - 0,9 V	0 [1/min]	0 [1/min]
1 V	5 [1/min]	24 [1/min]
10 V	120 [1/min]	600 [1/min]


Or control current 4 - 20 mAmps for rotational speed of 0 to max.

control current	PD 5101	PD 5106
0 - 4.9 mA	0 [1/min]	0 [1/min]
5 mA	5 [1/min]	24 [1/min]
20 mA	120 [1/min]	600 [1/min]

- 5 V TTL signal inverts sense of rotation (at standstill only !)
- A 5 volts TTL signal turns item ON, next signal to follow turns item OFF. This function is available all the time (in internal mode either).
- For triggering sense of rotation and ON / OFF function, a 5 volts pin is available in the interface connector, which allows potential-free switching.
- The analog interface is recognized at a level of 0.2 volts or a 4 m Amps current;

Interface socket wiring (37):

- Interface socket is located on the pump drive's rear panel; pin wiring:

PIN 1 = +5 V
PIN 2 = Start /Stop
PIN 3 = 0-10V
PIN 4 = sense of rotation
PIN 5 = 0 V(GND)
PIN 6 = 4-20 mA

Bild 47

view on terminal end

- Use IP 67 plug connectors only (6-pin 99-5121-15-06) Heidolph P/N 14-010-006-81 not to void splash protection.
- Make sure, plug connectors are installed properly.

9. PD 5201, PD 5206 (figs. 48, 49)

- PD 5201 and PD 5206 pump drives are splash-proof (IP 55).

Service work is limited to expert personnel authorized to do so by Heidolph.

Interface sockets (37 & 39) need to be closed with an interface plug or blind plug all the time, not to void splash protection.

9.1. Standard functions

- Turn item ON with power switch (29) located on the item's rear panel.
- Item runs through a processor self-test routine. The following parameters are displayed, one after the other:
 - all display (41) bars and all Select LEDs are ON.
 - pump head selected is displayed (e.g. „Pu 1”).
 - gearbox transmission ratio is displayed, e.g. „5.066”
 - hose ID selected is displayed, e.g. „3,1”
 - possible backstroke is displayed. e.g. „r on”.
- Item is ready to go at this point.
- Use sense of rotation key (31) to select sense of


fig. 48


Bild 49

- rotation (go by arrow marks). Sense of rotation may be changed on standstill only.

CAUTION: hit STOP key (32), select sense of rotation, START item

- Sense of rotation is confirmed by an LED (40).

- Start pump with „Start / Stop“ key (32). Pump running is confirmed by an LED (44).
- Hit „Start Stop“ key (32) once more to stop pump.
- When hitting the „max.“ key (33), pump will run at full speed as long as key (33) is held depressed. This feature speeds up filling and draining of hoses without doing any changes to the basic setting. "Full" as plain text appears in the display (41), when "max." key is depressed.
- Actual speed is displayed (41) in rpm (1/min), as long as the "speed rpm" LED is lighted.
- To alter speed setting (+ & -), use arrow keys (35).
- Feed rate is proportional with speed. Speed is displayed in rpm (1/min). For speed setting make reference to the feed rate graphs (Specifications).

9.2. Charging / Discharging of hose

As long as key "max" (33) is pressed pump runs with maximum speed. This function is useful to speed up the charging and discharging of the hose without changing the settings. When key "max" is pressed display shows "Full".

9.3. Extra functions

The PD 5201 and PD 5206 pump drives feature the following options:

- display of feed rate in ml/minute
- metering by volume
- interval metering
- calibrate feed rate or volume

To allow using the above options, you need to know about the pump model used; in case of peristaltic pumps, you need to know size of hose. In as shipped condition, the SP quick peristaltic (PU1) is preset. When using a different pump head, you need to program the item accordingly.

9.3.1. How to program particular pump heads (figs. 48,49)

If you have purchased a PD5201 or PD5206 in a product package, the drive has already been adjusted to the respective pump head at the factory!

- Turn power switch (29) OFF. Hit „Cal“ key (34) and turn power switch (29) ON. Pump head actually programmed is displayed.
- Select different pump head with arrow keys (35), e.g. „Pu 2“.
- Confirm choice with „Select“ key (36).
- For coding of pump heads refer to reference table below.

code	pump head	PD 5206	PD 5201
Pu 0	multi-channel pump head C 8		X
Pu 1	SP quick / SP quick d	X	X
Pu 2	SP Standard / SP vario	X	X
Pu 3	SP mini (obsolete)	X	X
Pu 4	Masterflex	X	X
Pu 8	multi-channel pump head C 12		X
Pu 9	multi-channel pump head C 4		X

9.3.2. Select display mode

Select key (36) allows to choose from several display modes:

- **speed rpm** (speed in 1/min)
- **tube mm** (hose ID in mm)
- **flow ml/min** (feed rate in ml/min)
- **vol ml** (metered volume in ml)
- **pause sec/min/h** (pause in seconds, minutes, or hours) (interval metering)

9.3.3. Display feed rate in ml/min

- In order to display feed rate in ml/min, the item needs to "know" about de hose diameter used.
- Feed rate is proportional with hose ID and pump speed.
- Enter hose ID: hit select-key (36) several times until „Tube mm“ LED lights.
- Enter hose ID with arrow keys (35) (e.g. 0.8; 1.7; 3.1; 4.8; 6.3 or 7.9).

If the LED **Vol ml** blinks when you adjust the internal diameter of the hose, the metering volume is being adapted to the minimum allowable volume. (See Chap. 9.3.5)

- Hit select key (36) until LED „Flow ml/min“ LED lights.'
- Feed rate in ml/min is displayed (41).
- To alter feed rate, use arrow keys (35) after having selected the "flow ml/min" mode with the select key.
- At this point, the processor is going to determine rotational speed.

- Since every hose has its own tolerance, it will be required to calibrate the pump drive for a particular hose, if high accuracy is mandatory.
- For geared pumps, the „tube mm“ menu is skipped, and „---“ displayed instead.

Caution: with reference to the very small quantities involved, feed rate is displayed in $\mu\text{l} / \text{min}$, when " **Cassette small** and hose size **0.2 0.5 or 0.9** are selected.

9.3.4. Calibrate feed rate

- Enter hose data and feed rate wanted, as described under 10.2.3.
- Hit „Cal“ key (34). Display (41) flashes.
- Hit „Start Stop“ key (32).
- The pump is going to run for 1 minute exactly. Collect fluid and determine quantity.
- Enter value so determined with arrow keys (35).
- Hold "Cal" key (34) depressed for not less than 3 seconds to save value. Data in memory are revised. "Cal" is displayed (41) for about 2 seconds to confirm action.
- Select actual feed rate wanted with arrow keys (35).

Caution: with reference to the very small quantities involved, feed rate is displayed in $\mu\text{l} / \text{min}$, when " **Cassette small** and hose size **0.2 0.5 or 0.9** are selected.

9.3.5. Metering by volume

There is a minimum allowable metering volume that depends upon the set speed and the hose diameter. Lower values can only be achieved with smaller hoses. (see page 178)

The PD5206 drive can only meter with a maximum of 300 l/min for reasons of precision! During the metering process, the drive is automatically limited to 300 l/min. The LED **Speed rpm** and the display blink before the metering process as an indication.

- The pump drive features metering of a certain volume of fluid. Once set, this quantity may be fed no matter how often by just hitting the Start/Stop key (32) or by triggering the interface (37/39). Max volume is 9999 ml
- To select a certain volume, you first need to enter hose ID.
- Feed rate is proportional to the number of rotations and hose ID.
- Enter hose ID: hit select-key (36) until „tube mm“ LED lights.
- Enter actual hose ID with arrow keys (35) (e.g. 0.8; 1.7; 3.1; 4.8; 6.3 or 7.9).

If the LED **Vol ml** blinks when you adjust the internal diameter of the hose, the metering volume is being adapted to the minimum allowable volume!

- Hold select-key (36) depressed until „Vol ml“ LED lights.

Caution: with reference to the very small quantities involved, feed rate is displayed in $\mu\text{l} / \text{min}$, when " **Cassette small** and hose size **0.2 0.5 or 0.9** are selected.

- Volume is displayed (41) in ml.
- Desired volume can be selected with arrow keys (35).
- At this point the processor is going to determine the number of revolutions required.
- Since every hose has its own tolerance, it will be required to calibrate the pump drive for a particular hose, if high accuracy is mandatory.

9.3.6. Calibrate metered volume

- Determine and set speed by your actual requirements

The PD5206 drive can only meter with a maximum of 300 l/min for reasons of precision! During the metering process, the drive is automatically limited to 300 l/min. The LED Speed rpm and the display blink before the metering process as an indication.

- Enter hose data and volume to be metered, as described under 10.3.5.
- Hit „Start Stop“ key (32).
- Pump is going to feed fluid at the speed and volume entered before.
- Determine exact quantity of fluid collected.
- Hit „Cal“ key (34). Display (41) flashes.

Caution: with reference to the very small quantities involved, feed rate is displayed in $\mu\text{l} / \text{min}$, when "Cassette small" and hose size **0.2 0.5 or 0.9** are selected.

- Enter value so determined with arrow keys (35).
- Hit "Cal" key (34) and hold depressed for not less than 3 seconds to enter. Data in memory are revised. "Cal" is displayed (41) for about 2 seconds to confirm action.
- Enter volume wanted with arrow keys (35).

Programmed flow rate and hence volume metered, are a function of pump head and hose selected. If, during calibration, the user finds a deviation in excess of 30 %, a setting error is very probable. In this case discontinue calibration (hit "select" key once), and go ahead checking the following parameters:

- did you select the right pump head (Chapter 10.3.1)?
- did you choose the right hose (Chapter 10.3.3)?
- is the hose properly installed in the pump head (Chapter 3) ? Is the clamping action ok ? (Chapter 3)

If wrong parameters have been entered for calibration of metered volume, metering action cannot be started with the "Start" key. The select LED will turn to „Flow ml/min“, and 4 dashes (----) are displayed.

Correct wrong parameters as described below:

- go into feed rate graph (Chapters 17 & 18) and determine a speed figure (1/min) plus a logical feed rate (ml/min).
- go to "speed" menu (1/min) with select key and choose a speed with arrow keys.
- go to "flow" menu (ml/min) with select key and choose a feed rate from the graph, as described under Chapter 10.3.4.
- at this point, a reasonable feed rate is saved.

9.3.7. Interval metering

- Interval metering is a sub-routine of metering by volume mode; repetition is not launched by hitting a key, or by a triggering pulse, but automatically after a pause selected before.
- Set volume as described under 10.3.5, metering by volume.
- Hold Select-key (36) depressed until „Pause sec/min/h" LED lights.
- Enter pause with arrow keys (35); choose pause between 0.1 seconds and 750 hours.
seconds range: display 0.1 to 599.9 (no extra signal displayed)
minutes range: display 10 to 599, „." displayed(600sec= 10min)
hours range: display 10 to 750 „h" displayed (600min=10Std.)
- Hold select key (36) depressed until „Vol ml" LED and LED „Pause sec/min/h" LED light the same time.
- Interval metering is controlled with the „Start Stop" key (32) or via the interface.
- After having hit the „Start Stop" key (32) the metered volume is displayed (41), rising from zero to the current value, or a pause interval rising from zero the value set before. The respective LED is going to flash.

9.3.8. Adjust backstroke for metering by volume mode (figs. 49)

To avoid dripping in the "metering by volume" mode, a short backstroke of the pump may be programmed.

- Turn power switch (29) OFF. Hit „Cal" key (34), and, while holding depressed, turn power switch (29) ON.
- Pump head, e.g. "Pu 1" is displayed (41).
- Hit arrow keys (35) until „r on" or „r off" is displayed (41).
- „r on" turns backstroke ON, „r off" turns backstroke OFF.
- Confirm mode with select key (36).

9.4. Interfaces (fig. 50)

9.4.1. Analog interface

- An analog interface connector (37) is installed in the item's rear panel.
- External pump control uses this analog interface.

Control voltage of 0 - 10 volts equals rotational speed between 0 and max.

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Control current 4 - 20 mA equals rotational speed between 0 and max.

control current	PD 5201	PD 5206
0 - 4.9 mA	0 [1/min]	0 [1/min]
5 mA	5 [1/min]	24 [1/min]
20 mA	120 [1/min]	606 [1/min]

- TTL-level 0 V = CW rotation, TTL-level 5 V = CCW rotation.
- A 5 volts TTL signal turns item ON, next signal to follow turns item OFF. This function is available all the time (in internal mode either).
- For triggering sense of rotation and ON / OFF function, a 5 volts pin is available in the interface connector, which allows potential-free switching.
- The analog interface is recognized at a level of 0.2 volts or a 4 m Amps current; This mode is confirmed by "EA" in the display. As long as the pump drive is controlled via the analog interface, Start / Stop-key will be the only active key on the item.
- Wiring of the interface connector (37):


fig. 49


PIN 1 = +5 volts
PIN 2 = start/stop
PIN 3 = 0-10 V
PIN 4 = sense of rotation
PIN 5 = 0 V (GND)
PIN 6 = 4-20 mA

View on terminal end

- Use IP 67 plug connectors only (Binder 6-pin 99-5121-15-06), Heidolph P/N 14-010-006-81, not to void splash water protection.
- Installation of plug connectors is limited to skilled personnel.

9.4.2. RS 232 interface (fig. 50)

- The RS 232 interface allows pump control from a PC via parallel interface.
- The RS 232 interface connector is installed in the rear panel of the pump drive. Wiring:

PIN 1 = PIN 1
PIN 2 = PIN 6
PIN 3 = PIN 2
PIN 4 = PIN 7
PIN 5 = PIN 3
PIN 6 = PIN 8
PIN 7 = PIN 4
PIN 8 = PIN 5


fig. 50

View on terminal (soldered) end

- Use IP 67 plug connectors only (Binder 8-pin 99-5171-15-08), Heidolph P/N 14-010-006-82, not to void splash water protection.
- Installation of plug connectors is limited to skilled personnel.

Interface parameters: 9600 Baud, no parity, 8 bit, 1 stop bit.

- As long as the RS 232 interface is active, all of the pump's control elements and displays will remain active.

Commands:

query	range	reply	remark	action
LEDx1?	x1=0..7	LED=0000 LED=0001	LED is „OFF" LED is „ON"	LED-status
DSP?		DSP=xxxx	xxxx=display	act. displays

command	range	action
TAx2!	x2=0..7 0=no key, 2=start, 4=max, 6=down, 7=select	key 1=Cal, 3=sense of rotation, 5=up, 7=select
SPx3=xxxx!	x3=0, xxxx=0001..5999(sec) x3=2, xxxx=0010..0599(min) x4=3, xxxx=0010..0750(hour)	pause
SDM=xxxx!	xxxx=000,1..9999	volume, metered
SDZ=xxxx!	xxxx=0000..9999 (depends on transmission used)	rotational speed
SMM=xxxx!	xxxx=000,1..9999(not during calibration) wait a minute for reply.	ml/min
SSD=xxxx!	xxxx=0000..0014	hose #

general	reply	remark	action
	OK	command successful	self explaining
	ERROR	command failed	none
	PO	metering complete	none

9.5. Recommended sizes of hose

- To achieve high metering accuracy, in particular in the metering by volume mode, we recommend the smallest hose diameter available; higher metering accuracy of small dia. hoses is a function of the higher number of revolutions of the pump drive.

C. RECOMMENDATIONS ABOUT HOSES

10. Hose properties

Designation	Tygon® Standard	Tygon® Food	Tygon® Hydrocarbons
Advantages	Low-price tube for laboratory applications, transparent. Ideal for inorganic media. Only slightly pervious to gas. Recommendable for viscous media. Good dielectric properties. Odourless, nontoxic.	Specially for food. Smooth inner walls, abrasion-proof. Neutral to smell and taste. Easy to clean and dry. Good dielectric properties.	Specially for hydrocarbons. Natural oil products and distillates. Ideal for petrol, kerosene, heating oil, cutting solutions, glycol-based coolants.
Limitations	Plasticizers may become disassociated. Limited life span.	Limited life span.	Unsuitable for concentrated acids and lyes. Plasticizers may become disassociated.
Use with: Acids Lyes Solvents Pressure Vacuum Viscous media Sterile media	good good unsuitable good good excellent under certain conditions	good good unsuitable good good excellent good	good good under certain conditions good good excellent under certain conditions
Temperature range	-50 to +70 °C	-40 to +70 °C	-40 to +70 °C
Physical properties	Thermoplastic soft PVC Transparent	Thermoplastic soft PVC Transparent	Thermoplastic soft PVC Yellow translucent
Complies with standards		FDA21 CFR 177.2600 BASPP and NSF (Standard 51) USDA Standard	
Permeability*)	CO₂ H ₂ O₂ N₂	1,7-7,4 97 0,7-12 0,2-3,0	1,7-7,4 97 0,7-12 0,2-3,0
Cleaning/sterilization	Can be sterilized with ethylene oxide or by autoclaving. Wrap the tube in a non-linting cloth or in paper and autoclave for 30 minutes at 121°C and 1 bar (takes on a milky colour). Dry in air at max. 66°C for 2 to 2 1/2 hours until the tube is clear again.	Can be sterilized with ethylene oxide or by autoclaving. Wrap the tube in a non-linting cloth or in paper and autoclave for 30 minutes at 121°C and 1 bar (takes on a milky colour). Dry in air at max. 66°C for 2 to 2 1/2 hours until the tube is clear again.	Sterilization not recommended.

Designation	PharMed®	Silikon	Viton®
Advantages	Ideal for non-toxic and non-haemolytic cell cultures and tissues. Impervious to normal light and UV radiation. The tube can be welded, bonded and formed. Very long life. Low gas permeability.	Excellent bio-compatibility. Practically no disassociation of plasticizers and additives. Neutral to smell and taste. Not susceptible to mildew. Ideal for low temperatures. Waterproof; resistant to ozone, radiation and sunlight. Does not warp. Silicon peroxide; hot vulcanized.	Highly resistant to chemicals. Excellent resistance to corrosive media, solvents and oils at high temperatures. Only slightly pervious to gas.
Limitations	Additives may become disassociated.	Unsuitable for concentrated acids, solvents, oils, acids or diluted caustic soda. Relatively highly pervious to gas.	Limited life span.
Use with: Acids Lyes Solvents Pressure Vacuum Viscous media Sterile media	good good unsuitable good excellent good excellent	under certain conditions under certain conditions unsuitable satisfactory good satisfactory excellent	excellent excellent varies, tests recommended good good good satisfactory
Temperature range	-60 to +130 °C	-50 to +230 °C	-30 to +200 °C
Physical properties	Thermoplastic elastomer on a polypropylene basis with plasticizers. Excellent tensile strength. Opaque, beige.	Polydimethylsiloxane with siliceous earth and silicon additives. Translucent, white. Excellent resistance to contact pressure.	Fluorocarbon rubber, thermally formed. Viton B (67 fluorinated). Opaque, black.
Complies with standards	USP, Class VI FDA21 CFR 177.2600 NSF (standard 51)	USP, Class VI FDA21 CFR 177.2600	
Permeability*) CO₂ H₂ O₂ N₂	450 - 75 29	20132 6579 7961 2763	76-79 - 13-15 4,3
Cleaning/sterilization	Can be autoclaved without ageing symptoms.	Can be cleaned with hot water and suds. Use only soap without oil additives; use no synthetic detergents. Rinse with distilled water. Can be sterilized by autoclaving in a moist atmosphere or with gamma rays. Cannot be sterilized with ethylene oxide.	Can be sterilized at 249°C in 16 hours in a circulating air heating cabinet.

11. Chemical stability of hoses

- Tube:

P = PharMed®
 S = Silikon
 T = Tygon® standard und Tygon® food
 K = Tygon® hydrocarbon
 V = Viton®

- Resistance

A = very good
 B = good
 C = satisfactory
 D = unsuitable
 - = not tested

- Please note:


- All information provided is not guaranteed
- Recommend testing of tubing prior to application use

Chemical	P	S	T	K	V	Chemical	P	S	T	K	V
A Acetaldehyde	D	C	D	D	D	Bromine	D	D	D	D	A
Acetic acid, 10% in W.	A	A	A	A	-	Butane	A	A	A	A	A
Acetic acid, 100%	B	D	D	D	-	Butanol (butyl alcohol)	D	D	D	A	A
Acetic anhydride	A	A	D	D	D	Butyl acetate	B	D	D	D	D
Acetone	D	C	D	D	D	Butyric acid	B	D	D	C	-
Acetonitrile	D	D	D	D	D	C Calcium oxide	A	A	A	A	-
Acetyl bromide	C	D	D	D	-	Carbon bisulfide	D	D	D	D	-
Acetyl chloride	C	D	D	D	A	Carbon tetrachloride	D	D	D	A	A
Aliphatic hydrocarbons	D	D	D	B	-	Chlorine, wet	D	D	B	B	B
Aluminium chloride, 53% in W.	A	A	A	A	A	Chloracetic acid, 20% in W.	B	A	A	D	D
Aluminium sulfate, 50% in W.	A	A	A	A	A	Chlorobenzene	D	D	D	D	A
Alums	A	A	A	A	-	Chrobromomethane	B	D	D	D	A
Ammonia, gas and liquid	A	D	B	B	D	Choroform	D	D	D	D	A
Ammonium acetate, 45% in W.	A	A	A	A	-	Chromic acid, 20% in W.	A	D	B	C	A
Ammonium carbonate, 20% in W.	A	A	A	A	-	Chromic acid, 50% in W.	C	D	C	D	-
Ammonium chloride	A	C	A	A	A	Copper salts	A	A	A	A	-
Ammonium hydroxide, 30% in W.	A	D	A	C	A	Cyclohexane	D	D	D	C	A
Ammonium nitrate	A	C	A	A	-	Cyclohexanone	D	D	D	D	D
Ammonium phosphate	A	A	A	A	-	Cyrosulfonic acid	D	D	D	D	D
Ammonium sulfate	A	A	A	A	A	D Diesel	D	D	-	B	-
Amyl acetate	C	D	D	D	D	Dimethyl formamide	B	A	D	D	D
Amyl alcohol	D	D	D	A	A	E Ethanol (ethyl alcohol)	A	B	D	B	A
Amyl chloride	C	D	D	D	-	Ether	C	D	D	C	-
Aniline	C	D	D	D	-	Ethyl acetate	B	D	D	D	D
Aniline hydrochloride	C	D	D	D	D	Ethyl bromide	D	D	D	D	-
Aqua regia 80% HCl, 20% HNO ₃	D	D	D	D	-	Ethyl chloride	C	D	D	D	A
Aromatic hydrocarbons	D	D	D	D	-	Ethylamine	D	C	D	D	-
Arsenic salts	A	A	A	A	-	Ethylene chlorhydrin	A	B	D	B	A
B Barium salts	A	A	A	A	-	Ethylene dichloride	C	D	D	D	B
Benzaldehyde	D	C	D	D	D	Ethylene glycol	A	A	A	A	A
Benzene	D	D	D	D	-	Ethylene oxide	A	A	A	A	D
Benzenesulfonic acid	D	D	D	D	A	F Ferric chloride 40% in W.	A	A	A	A	B
Boric acid, 4% in W.	A	A	A	A	A	Ferric sulfate 5% in W.	A	A	A	A	A

Chemical	P	S	T	K	V	Chemical	P	S	T	K	V
Ferrous chloride 43% in W.	A	A	A	A	-	Oleic acid	C	B	D	B	B
Ferrous sulfate 5% in W.	A	A	A	A	-	Perchloric acid. 67% in W.	A	D	C	D	A
Fluoboric acid, 10% in W.	D	D	A	A	-	Perchlorethylene	C	D	D	D	A
Fluoroborate salts	A	-	A	A	-	Phenol, 91% in W.	A	B	D	C	-
Fluosilicicacid	C	B	D	B	-	Phosphonic acid, 25% in W.	A	D	A	A	A
Formaldehyde, 37% in W.	D	C	D	D	D	Phthalic acid, 9% in Alc.	A	B	D	C	-
Formic acid, 25% in W.	A	A	A	C	D	Potassium carbonate, 55% in W.	A	A	A	A	-
Freon 11	A	A	A	A	-	Potassium cyanide, 33% in W.	A	A	A	A	-
Fruit juice	A	A	A	A	A	Potassium hydroxide,<10% in W	A	A	A	D	B
G Gasoline, high-aromatic	D	D	D	B	A	Potassium iodide, 56% in W.	A	A	A	A	-
Gasoline, non-aromatic	D	D	D	B	A	Propanol (propyl alcohol)	C	A	D	D	B
Glycerin	A	A	A	A	A	Pyridine	C	D	D	D	D
H Hydrobromic acid, 20-50%	A	A	A	A	C	S Silicone oils	C	D	B	A	A
Hydrochloric acid, 10% in W.	A	A	A	A	A	Silver nitrate, 55% in W.	A	A	A	A	A
Hydrochloric acid, 37% in W.	B	B	A	A	-	Soap solutions	B	A	A	A	A
Hydrocyanic acid	A	A	A	A	A	Sodium bicarbonate, 7% in W.	A	A	A	A	A
Hydrofluoric acid, 10% in W.	D	D	A	A	B	Sodium bisulfate	A	-	A	A	-
Hydrofluoric acid, 50%	D	D	B	D	A	Sodium borate	A	A	A	A	A
Hydrogen peroxide, 10% in W.	A	A	A	A	A	Sodium carbonate	A	A	A	A	B
Hydrogen peroxide, 90% in W.	B	C	D	D	-	Sodium ferrocyanide	A	-	A	A	-
Hydroiodic acid	B	B	A	A	-	Sodium hydrosulfite	A	-	A	A	-
Hypochlorousacid,25% in W.	A	A	A	A	A	Sodium hydroxide, 10-15% in W.	A	A	A	D	B
I Iodine solutions	A	C	A	A	-	Sodium hydroxide, 30-40% in W.	A	C	A	D	B
K Ketones	D	D	D	D	-	Sodium nitrate, 3.5% in W.	A	A	A	A	-
L Lactic acid, 10% in W.	A	A	A	A	-	Sodium sulfate, 3.6% in W.	A	A	A	A	A
Lactic acid, 85% in W.	B	D	D	D	-	Sodium sulfide, 13% in W.	A	A	A	A	-
Lead acetate, 35% in W.	A	A	A	A	-	Stearic acid, 5% in Alc.	C	D	D	B	-
M Maqanese salts	A	A	A	A	-	Sulfuric acid, 10% in W.	A	A	A	B	A
Magnesium chloride, 35% in W.	A	A	A	A	A	Sulfuric acid, 30% in W.	A	B	A	B	A
Magnesium sulfate, 25% in W.	A	A	A	A	-	Sulfuric acid, 95-98% in W.	D	D	D	A	
Mercury salts	A	A	A	A	-	Sulfurous acid	A	A	A	A	A
Methane	A	-	A	A	A	T Tannic acid, 75% in W.	B	A	B	D	-
Methanol	A	B	D	B	D	Tartaric acid, 56% in W.	A	A	A	A	A
Methyl Ethyl Ketone	D	D	D	D	D	Tin salts	A	A	A	A	-
Monoethanolamine	C	D	D	D	D	Toluene (toluol)	D	D	D	D	A
N Naphtha	D	D	D	D	A	Trichloroacetic acid, 90% in W.	B	A	A	D	C
Nickel salts	A	A	A	A	-	Trichlorethylene	C	D	D	D	A
Nitric acid, 10% in W.	A	C	A	D	A	Trisodium phosphate	A	A	A	A	A
Nitric acid, 35% in W.	A	D	A	D	A	Turpentine	D	D	D	B	A
Nitric acid, 68-71% in W.	D	D	D	D	-	U Urea, 20% in W.	A	A	A	A	-
Nitrobenzene	D	D	D	D	-	Uric acid	A	A	A	C	-
Nitrous acid, 10% in W.	A	B	A	C	-	X Xylene	D	D	D	D	B
O Oil, animal	C	A	D	A	-	Z Zinc chloride, 80% in W.	A	A	A	A	A
Oils, mineral	D	D	C	A	A						

12. Service life of hoses

Peristaltic pumphoses are wear parts. Service life of hoses is a function of material, rotational speed, and properties of the substance to be fed, e.g. temperature, viscosity, chemical composition, and pressure. For more information about expected service life of hoses (not binding), refer to graphs overleaf.

13. Feed rates - list of hoses for peristaltic pumps

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14. Feed rates - list of hoses for multi-channel pumps

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Flow rates in combination with pump head / pump drive:

PD 5101 / PD 5201		min.	max.	min.	max.	min.	max.	min.	max.	max. number	
Cassette medium / Pump head C8	[ml/min]	0,29	7,20	1,0	25,2					8	
Cassette large / Pump head C8	[ml/min]			1,0	25,2	3,8	96,0	8,2	204	13,5 336	
PD 5001											
Cassette medium / Pump head C8	[ml/min]	0,60	7,20	2,10	25,2					8	
Cassette large / Pump head C8	[ml/min]			2,10	25,2	8,00	96,0	17,0	204	28,0 336	

Tubing order numbers:

Silikon	Meterware	525-33000-00	525-34000-00	525-36000-00	525-30027-00	525-30028-00
Viton	Meterware	525-53000-00	525-54000-00	525-56000-00	525-50027-00	525-50028-00
PharMed	Meterware	525-23000-00	525-24000-00	525-26000-00	525-20027-00	525-20028-00
Tygon (Standard)	Meterware	525-63000-00	525-64000-00	525-66000-00	525-60027-00	525-60028-00
Tygon (Hydrocarbons)	Meterware	525-73000-00	525-74000-00	525-76000-00	525-70027-00	525-70028-00
Tygon (Food)	Meterware	525-83000-00	525-84000-00	525-86000-00	525-80027-00	525-80028-00

Cassette / Pump head / Pump drive order numbers:

Pump drive PD 5001 (incl. Adaptor)	523-50013-00	Multi-channel pump head C12	524-81220-00	Cassette small	524-90022-00
Pump drive PD 5101 (incl. Adaptor)	523-51013-00	Multi-channel pump head C4	524-80420-00	Cassette medium	524-90021-00
Pump drive PD 5201 (incl. Adaptor)	523-52013-00	Multi-channel pump head C8	524-40810-00	Cassette large	524-90010-00

CLEANING AND SERVICING

For **cleaning** wipe housing surfaces with a damp cloth (add some sort of mild liquid soap).

NOTE: To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning pads containing metallic particles.

Pump drives and pump heads are maintenance-free. Repair work is limited to technicians so authorized by Heidolph. Please call your local Heidolph dealer or field representative.

Hoses for peristaltic pumps are wear parts. We recommend to replace hoses in regular intervals to avoid the mess with broken hoses. Service life of hoses used in conventional peristaltic pumps is between about 50 and 500 hours in use (depends on hose material and mission profile).

UNINSTALL, FORWARD & STORE

Uninstall

1. Turn the item OFF and disconnect mains plug.
2. Clear area around the item to get enough room for uninstall work.
3. Uninstall optional equipment.

Forward & Store

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

Caution
Do not jolt or shake the item during transport.

DISPOSAL

For disposal, please comply with your local or national regulations.

Split by metal, plastic, etc.

Packing material to be treated as described above (material split).

TROUBLESHOOTING

Work on electric and electronic components is limited to qualified personnel.

pump won't feed, motor running

- a) SP quick: hose bed open
- b) hose with deviating wall thickness installed (refer to Chapter 3.1)
- c) hose with deviating wall thickness tolerance installed; use genuine Heidolph hose only.
- d) SP vario: wall thickness of hose not set (refer to Chapter 3.6)
- e) pump head (SP Standard, SP mini or SP vario) not properly seated on drive, coupling not in gear

Hose creeping in peristaltic pump

- a) re-adjust hose clamp (refer to Chapter 3.3.)
- b) use short hoses for 0.8 mm and 1.7 mm hoses for clamping (refer to Chapter 3.7).

Pump drive won't run

- a) wall outlet dead
- b) check power cord and pump drive wiring
- c) fuse blown (replace fuse in rear panel of PD 5001/5006 (use 1.25 Amps delay-type fuse))
- d) power switch broken
- e) control circuit or drive motor broken


SPECIFICATIONS

15. Peristaltic pump drives


model	PD 5001	PD 5006	PD 5101	PD 5106	PD 5201	PD 5206				
P/N	523-50010-00	523-50060-00	523-51010-00	523-51060-00	523-52010-00	523-52060-00				
electronic motor control	capacitor motor	Electronically commutated direct current (DC) motor								
feed rates	refer to Chapters 16 & 17									
speed range (1/min.)	10-120	50-600	5-120	24-600	5-120	24-600				
speed control	rotary knob, scaled				4-digit digital display					
speed governor	analog with answer-back		digital, with answer-back							
accuracy of control	±2%		± 0,5%							
sense of rotation	CW, CCW selector									
motor power (W)	71		100							
dissipated power (W)	150		140							
interface, analog			speed: (0-10 V or 4-20 mA), sense of rotation & Start/Stop (TTL)							
interface, digital					RS232					
feed rate display					4-digit digital display					
metering by volume (ml)					0.1-9999					
interval metering (ml)					0.1-9999; pause 0., sec-750 h					
motor protector	ÜHS*)	ÜHS*)	Electronic current limitation and overheat protection (ÜHS*))		Electronic current limitation and overheat protection (ÜHS*))					
weight, pump head not included (kg)	6.8	6.2	8.2	7.6	8.3	7.7				
protective class	IP30		IP55							
ambient conditions	0-40 °C at 80 % rel. humidity, no dewing									
size (w x d x h) (mm)	166x256x225									

*) ÜHS = overheat protection

16. Feed rate graphs for peristaltic pumps

17. Feed rate graphs for multi-channel pumps

WARRANTY, LIABILITY & COPYRIGHT

Warranty

Heidolph Instruments provides a three-year warranty on the products described here (with the exception of glass and consumable parts) if registered with enclosed warranty card or via internet (www.heidolph.com). Warranty starts with the date of registration. Without registration warranty starts according to serial number. This warranty covers defects in materials and workmanship. Damage in transit is excluded from this warranty.

To file for such warranty service, contact Heidolph Instruments (phone +49-9122-9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty.

Alterations to the present warranty need Heidolph Instruments' consent in writing.

Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

FAQ / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For repair services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

Note

You will receive approval for sending your defective item to the following address:

Heidolph Instruments GmbH & Co. KG
Lab Equipment Sales
Walpersdorfer Str. 12
D-91126 Schwabach / Germany
phone: ++49-9122-9920-68
Fax: ++49-9122-9920-65
E-Mail: Sales@Heidolph.de

Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take proper *protective* meason to ensure the safety of our receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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P S T K V

A	Acetaldehyde	D	C	D	D
	Acetic acid, 10% in W.	A	A	A	-
	Acetic acid, 100%	B	D	D	-
	Acetic anhydride	A	A	D	D
	Acetone	D	C	D	D
	Acetonitrile	D	D	D	D
	Acetyl bromide	C	D	D	-
	Acetyl chloride	C	D	D	A
	Aliphatic hydrocarbons	D	D	B	-
	Aluminium chloride, 53% in W.	A	A	A	A
	Aluminium sulfate, 50% in W.	A	A	A	A
	Alums	A	A	A	-
	Ammonia, gas and liquid	A	D	B	D
	Ammonium acetate, 45% in W.	A	A	A	-
	Ammonium carbonate, 20% in W.	A	A	A	-
	Ammonium chloride	A	C	A	A
	Ammonium hydroxide, 30% in W.	A	D	A	C
	Ammonium nitrate	A	C	A	-
	Ammonium phosphate	A	A	A	-
	Ammonium sulfate	A	A	A	A
	Amyl acetate	C	D	D	D
	Amyl alcohol	D	D	A	A
	Amyl chloride	C	D	D	-
	Aniline	C	D	D	-
	Aniline hydrochloride	C	D	D	D
	Aqua regia 80% HCl, 20% HNO ₃	D	D	D	-
	Aromatic hydrocarbons	D	D	D	-
	Arsenic salts	A	A	A	-
B	Barium salts	A	A	A	-
	Benzaldehyde	D	C	D	D
	Benzene	D	D	D	-
	Benzenesulfonic acid	D	D	D	A
	Boric acid, 4% in W.	A	A	A	A
	Bromine	D	D	D	A

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P S T K V

Butane	A	A	A	A	
Butanol (butyl alcohol)	D	D	D	A	
Butyl acetate	B	D	D	D	
Butyric acid	B	D	D	C	
C	Calcium oxide	A	A	A	-
Carbon bisulfide	D	D	D	-	
Carbon tetrachloride	D	D	D	A	
Chlorine, wet	D	D	B	B	
Chloracetic acid, 20% in W.	B	A	A	D	
Chlorobenzene	D	D	D	A	
Chrobromomethane	B	D	D	A	
Choroform	D	D	D	A	
Chromic acid, 20% in W.	A	D	B	C	
Chromic acid, 50% in W.	C	D	C	D	
Copper salts	A	A	A	A	
Cyclohexane	D	D	C	A	
Cyclohexanone	D	D	D	D	
Cyrosulfonic acid	D	D	D	D	
D	Diesel	D	D	-	B
Dimethyl formamide	B	A	D	D	
E	Ethanol (ethyl alcohol)	A	B	D	B
Ether	C	D	D	C	
Ethyl acetate	B	D	D	D	
Ethyl bromide	D	D	D	-	
Ethyl chloride	C	D	D	A	
Ethylamine	D	C	D	-	
Ethylene chlorhydrin	A	B	D	B	
Ethylene dichloride	C	D	D	B	
Ethylene glycol	A	A	A	A	
Ethylene oxide	A	A	A	D	
F	Ferric chloride 40% in W.	A	A	A	B
Ferric sulfate 5% in W.	A	A	A	A	
Ferrous chloride 43% in W.	A	A	A	-	
Ferrous sulfate 5% in W.	A	A	A	-	

Produits chimiques	P S T K V	Produits chimiques	P S T K V
Fluoboric acid, 10% in W.	D D A A -	P Perchloric acid, 67% in W.	A D C D A
Fluoroborate salts	A - A A -	Perchlorethylene	C D D D A
Fluosilicicacid	C B D B -	Phenol, 91% in W.	A B D C -
Formaldehyde, 37% in W.	D C D D D	Phosphonic acid, 25% in W.	A D A A A
Formic acid, 25% in W.	A A A C D	Phthalic acid, 9% in Alc.	A B D C -
Freon 11	A A A A -	Potassium carbonate, 55% in W.	A A A A -
Fruit juice	A A A A A	Potassium cyanide, 33% in W.	A A A A -
G Gasoline, high-aromatic	D D D B A	Potassium hydroxide,<10% in W	A A A D B
Gasoline, non-aromatic	D D D B A	Potassium iodide, 56% in W.	A A A A -
Glycerin	A A A A A	Propanol (propyl alcohol)	C A D D B
H Hydrobromic acid, 20-50%	A A A A C	Pyridine	C D D D D
Hydrochloric acid, 10% in W.	A A A A A	S Silicone oils	C D B A A
Hydrochloric acid, 37% in W.	B B A A -	Silver nitrate, 55% in W.	A A A A A
Hydrocyanic acid	A A A A A	Soap solutions	B A A A A
Hydrofluoric acid, 10% in W.	D D A A B	Sodium bicarbonate, 7% in W.	A A A A A
Hydrofluoric acid, 50%	D D B D A	Sodium bisulfate	A - A A -
Hydrogen peroxide, 10% in W.	A A A A A	Sodium borate	A A A A A
Hydrogen peroxide, 90% in W.	B C D D -	Sodium carbonate	A A A A B
Hydroiodic acid	B B A A -	Sodium ferrocyanide	A - A A -
Hypochlorousacid,25% in W.	A A A A A	Sodium hydrosulfite	A - A A -
I Iodine solutions	A C A A -	Sodium hydroxide, 10-15% in W.	A A A D B
K Ketones	D D D D -	Sodium hydroxide, 30-40% in W.	A C A D B
L Lactic acid, 10% in W.	A A A A -	Sodium nitrate, 3.5% in W.	A A A A -
Lactic acid, 85% in W.	B D D D -	Sodium sulfate, 3.6% in W.	A A A A A
Lead acetate, 35% in W.	A A A A -	Sodium sulfide, 13% in W.	A A A A -
M Maqanese salts	A A A A -	Stearic acid, 5% in Alc.	C D D B -
Magnesium chloride, 35% in W.	A A A A A	Sulfuric acid, 10% in W.	A A A B A
Magnesium sulfate, 25% in W.	A A A A -	Sulfuric acid, 30% in W.	A B A B A
Mercury salts	A A A A -	Sulfuric acid, 95-98% in W.	D D D D A
Methane	A - A A A	Sulfurous acid	A A A A A
Methanol	A B D B D	T Tannic acid, 75% in W.	B A B D -
Methyl Ethyl Ketone	D D D D D	Tartaric acid, 56% in W.	A A A A A
Monoethanolamine	C D D D D	Tin salts	A A A A -
N Naphtha	D D D D A	Toluene (toluol)	D D D D A
Nickel salts	A A A A -	Trichloroacetic acid, 90 in W.	B A A D C
Nitric acid, 10% in W.	A C A D A	Trichlorethylene	C D D D A
Nitric acid, 35% in W.	A D A D A	Trisodium phosphate	A A A A A
Nitric acid, 68-71% in W.	D D D D -	Turpentine	D D D B A
Nitrobenzene	D D D D -	U Urea, 20% in W.	A A A A -
Nitrous acid, 10% in W.	A B A C -	Uric acid	A A A C -
O Oil, animal	C A D A -	X Xylene	D D D D B
Oils, mineral	D D C A A		
Oleic acid	C B D B B		

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